

Effect of bottle height on the corneal endothelium during phacoemulsification

Hisaharu Suzuki, MD, PhD, Kotaro Oki, MD, PhD, Toshihiko Shiwa, MD, PhD,
Hideaki Oharazawa, MD, PhD, Hiroshi Takahashi, MD, PhD

PURPOSE: To directly measure intraocular pressure (IOP) in simulated phacoemulsification and to assess the usefulness of lowering the bottle height in protecting the corneal endothelium in clinical phacoemulsification.

SETTING: Nippon Medical School Hospital, Tokyo, Japan.

METHODS: Simulated phacoemulsification was performed in porcine eyes with 2 bottle heights, 65.0 cm (BH 65 group) and 19.0 cm (BH 19 group). The IOP was continuously measured with a microprobe. In a clinical study, phacoemulsification was performed with a bottle height of 60.0 cm (BH 60 group) and of 30.0 cm (BH 30 group). One day, 1 week, and 1 and 3 months after surgery, cell density and corneal volumes were measured using specular microscopy and rotating Scheimpflug photography, respectively.

RESULTS: In the simulation study, IOP fluctuated between 50 mm Hg and 60 mm Hg in the BH 65 group and between 20 mm Hg and 30 mm Hg in the BH 19 group. In the clinical study of 31 eyes, the rate of cell density decrease was significantly lower in the BH 30 group than in the BH 60 group at all time points. The rate of increase in corneal volume was significantly lower in the BH 30 group than in the BH 60 group at 1 month.

CONCLUSIONS: Intraoperative IOP in phacoemulsification with a usual bottle height appeared to exceed the normal range. Phacoemulsification with a low bottle height was less harmful to the corneal endothelium.

J Cataract Refract Surg 2009; 35:2014–2017 © 2009 ASCRS and ESCRS

The corneal endothelium can be damaged during phacoemulsification by factors such as excessive duration of phacoemulsification^{1–3}; localized temperature increases⁴; or damage from lens nucleus fragments caused by the turbulent flow of the irrigating solution,^{5,6} air bubbles,^{7,8} or free radicals associated with ultrasound oscillation.^{9–13} However, little attention

has been paid to intraocular pressure (IOP) during surgery as a cause of corneal endothelial damage. Intraoperative IOP is determined by a balance between the irrigating (influx) pressure, which depends on the height of the irrigating bottle, and the outflow volume, which depends on the vacuum pressure, aspiration rate, and fluid leakage from the incision site. Thus, the intraoperative IOP at one moment can be considered an instantaneous measure of the fluid flow into the anterior chamber.

Although turbulent flow of the irrigating solution is a known cause of corneal endothelial damage,^{5,6} it has been discussed mainly in terms of the lens nucleus colliding with the corneal endothelium. The effect of elevated IOP or alteration of intraoperative IOP on the corneal endothelium has not been fully elucidated. In this study, we directly measured intraoperative IOP in simulated phacoemulsification in animal eyes and evaluated the effect of lowering the bottle height on corneal endothelial damage in clinical phacoemulsification.

Submitted: March 23, 2009.

Final revision submitted: May 15, 2009.

Accepted: May 17, 2009.

From the Department of Ophthalmology, Nippon Medical School, Tokyo, Japan.

No author has a financial or proprietary interest in any material or method mentioned.

Corresponding author: Hiroshi Takahashi, MD, PhD, Department of Ophthalmology, Nippon Medical School, 1-1-5 Sendagi, Bunkyo-ku, Tokyo 113-8602, Japan. E-mail: tash@nms.ac.jp.

MATERIALS AND METHODS

Simulation Study

Phacoemulsification was performed in porcine eyes obtained from a local abattoir using a Sovereign WhiteStar system (Abbott Medical Optics). The surgery was performed via a superior 2.4 mm sclerocorneal incision with an aspiration rate of 22 mL/min and a vacuum pressure of 150 mm Hg. Two bottle height settings, 65.0 cm (BH 65 group) and 19.0 cm (BH 19 group), were used. Because the bottle height value that is set and displayed by the phaco device is not the true distance from the eye but is the distance from a certain point of the device to the bottle, actual distances from the eye to the bottle (true bottle height) were measured and found to be 88.0 cm in the BH 65 group and 40.5 cm in the BH 19 group. Intraoperative IOP was measured using a universal pressure meter (DPM-2 Plus, Fluke Corp.). The microprobe of the device was inserted from the 6 o'clock limbus into the anterior chamber, and the IOP was recorded continuously during the procedure.

Clinical Study

The clinical study included patients who had phacoemulsification and intraocular lens implantation over a 3-month period starting in November 2005 at Nippon Medical School Hospital. In all cases, lens nucleus hardness was approximately 2 according to the Emery-Little classification.

The same surgeon (H.S.) operated on all patients using the same type of phaco system used in the simulation study. Phacoemulsification was performed through a superior sclerocorneal incision using sodium hyaluronate 2.3% (Healon5), an ultrasound power output of 20%, and vacuum pressure of 150 mm Hg. Two bottle height settings, 60.0 cm (BH 60 group) and 30.0 cm (BH 30 group), were used; both bottle heights were displayed on the device. These settings were different from those in the simulation study. The reason is that a bottle height of 19.0 cm was considered too low and there was concern that it may cause anterior chamber collapse in clinical phacoemulsification; in addition, the 60.0 cm height is the standard setting used at the hospital where the surgery was performed. The aspiration rate was set at 20 mL/min in the BH 60 group and 18 mL/min in the BH 30 group to prevent anterior chamber collapse.




Figure 1. Intraoperative IOP in the BH 65 and BH 19 groups. In the BH 65 group, intraoperative IOP oscillated between 50 mm Hg and 60 mm Hg and in the BH 19 group, between 20 mm Hg and 30 mm Hg (BH19 = bottle height at 19.0 cm; BH65 = bottle height at 65.0 cm; IOP = intraocular pressure).

A specular microscope (Noncon Robo, Konan Medical) and a rotating Scheimpflug camera (Pentacam, Oculus) were used to analyze the cornea before surgery and 1 day, 1 week, and 1 and 3 months after surgery. Specular microscopy was performed in the central cornea to measure cell density and to analyze polymegathism (coefficient of variation) and pleomorphism (hexagonal cell ratio 6A) using the center method. The Scheimpflug camera was used to calculate the corneal volume within 3.0 and 10.0 mm circles in the central cornea.^{14,15}

Statistical analysis was performed using the unpaired Student *t* test, F test, 1-way analysis of variance (ANOVA), or Dunnett test, when appropriate.

RESULTS

Simulation Study

The intraoperative IOP fluctuated in an oscillating fashion (Figure 1). The IOP in the BH 65 group (mean 52.72 ± 7.94 mm Hg) was statistically significantly higher than in the BH 19 group (mean 21.96 ± 4.98 mm Hg) ($P < .05$, Student *t* test). Moreover, the range between the maximum IOP and minimum IOP was significantly greater in the BH 65 group than in the BH 19 group ($P < .05$, F test).

Clinical Study

The clinical study comprised 63 patients (63 eyes); there were 31 patients in the BH 60 group and 32 patients in the BH 30 group. In all cases, the surgery lasted fewer than 20 minutes and the effective phaco time was generally fewer than 10 seconds.

At all time points, the rate of cell density decrease was statistically significantly lower in the BH 30 group than in the BH 60 group ($P < .05$, ANOVA and Dunnett test) (Figure 2). No significant difference was observed in the coefficient of variation and hexagonal cell ratio between 2 groups at any time point.

Figures 3 and 4 show the rate of corneal volume increase compared with the preoperative value. One




Figure 2. Rate of cell density decrease (* = $P < .05$, ANOVA and Dunnett test; BH30 = bottle height at 30.0 cm; BH60 = bottle height at 60.0 cm).




Figure 3. Rate of 3.0 mm corneal volume increase (* = $P < .05$, ANOVA and Dunnett test; BH30 = bottle height at 30.0 cm; BH60 = bottle height at 60.0 cm).

month after surgery, the corneal volumes within the 3.0 mm and 10.0 mm circles were significantly smaller in the BH 30 group than in the BH 60 group ($P < .05$, ANOVA and Dunnett test).

DISCUSSION

When the caliber or the resistance of the irrigating tube is not taken into consideration, the irrigating pressure for the eye can be roughly estimated by the following formula:

$$\text{Irrigating pressure} = \text{TBH}(\text{cm}) \times 10 / 13.5 (\approx 0.74) \text{ mmHg}$$

where TBH is the true bottle height and 13.5 is the specific gravity of mercury. Using animal eyes, we directly measured IOP during simulated phacoemulsification under 2 bottle heights, 65.0 cm (BH 65 group) and 19.0 cm (BH 19 group). The true bottle height was 88.0 cm in the BH 65 group and 40.5 cm in the BH 19 group. Thus, the rough estimate of irrigating pressure is 65 (88×0.74) mm Hg in the BH 65 group and 30 (40.5×0.74) mm Hg in the BH 19 group. Results of measurements using a microprobe for the most part agreed with these estimates and clearly showed that intraoperative IOP with the usual bottle height setting (as in the BH 65 group) was higher than the normal range. We also found that the intraoperative IOP fluctuated in a sawtooth wave pattern and that the range between the peak and trough of the wave was significantly greater in the BH 65 group than in the BH 19 group, suggesting that the higher the irrigating pressure, the more significant the alteration in the fluid flow in the anterior chamber.

In our clinical study, the postoperative corneal endothelial damage was significantly less with a low bottle height than with the usual height setting, as shown by parameters such as the rate of cell density increase and the increase in corneal volume, which we have proposed as a new parameter of corneal




Figure 4. Rate of 10.0 mm corneal volume increase (* = $P < .05$, ANOVA and Dunnett test; BH30 = bottle height at 30.0 cm; BH60 = bottle height at 60.0 cm).

endothelial damage after phacoemulsification.^{14,15} Unlike the cell density value, the change in corneal volume, with the 10.0 mm volume in particular, could reflect the total histologic alteration in the cornea, induced not only by endothelial damage but also by mechanical trauma at an incision site. We found little alteration in corneal volume with a lower bottle height, suggesting minimum damage to the entire cornea.

It is well known that in acute glaucoma, extremely high IOP can be hazardous to the corneal endothelium and that the duration of acute glaucoma episodes may be a risk factor for corneal endothelial damage.¹⁶ In our study, however, all surgery was completed within 20 minutes, suggesting that the effect of elevated IOP on the corneal endothelium during the surgery may not be significant. Results in the simulation study suggest that the beneficial effects of a lower bottle height is the result of a smaller fluctuation in intraoperative IOP, which implies a more subtle alteration in irrigating fluid flow than that caused by the usual bottle height. In addition, improved retention of the ophthalmic viscosurgical device (OVD) in the anterior chamber during surgery with a lower bottle height likely plays a role in protecting the corneal endothelium. In addition, with slow fluid flow, movements of the phacoemulsification probe are restrained. Probe movements can also contribute to damage at the incision site, as suggested by the corneal volume analysis. Furthermore, we lowered the aspiration rate in the BH 30 group to prevent anterior chamber collapse, which contributed to the slow fluid flow associated with lower irrigating pressure. In any case, our results clearly show that lowering the bottle height helps to protect the corneal endothelium. However, whether lowering intraoperative IOP contributes to the effect remains to be elucidated. To evaluate this, we are planning a histologic analysis to determine how intraoperative IOP affects the corneal endothelium.

In conclusion, phacoemulsification with a low bottle height was less harmful to the corneal endothelium than phacoemulsification with a normal bottle height. The beneficial effects of lowering bottle height may result from minimizing the impact of lens nucleus fragments with the corneal endothelium in the setting of milder turbulence in the anterior chamber, restrained movement of the phacoemulsification probe, better retention of the OVD in the anterior chamber, or decreased IOP.

REFERENCES

1. Olson LE, Marshall J, Rice NSC, Andrews R. Effects of ultrasound on the corneal endothelium: I. The acute lesion. *Br J Ophthalmol* 1978; 62:134–144. Available at: <http://www.pubmedcentral.nih.gov/picrender.fcgi?artid=1043167&blobtype=pdf>. Accessed July 14, 2009
2. Beesley RD, Olson RJ, Brady SE. The effects of prolonged phacoemulsification time on the corneal endothelium. *Ann Ophthalmol* 1986; 18:216–219; 222
3. Pirazzoli G, D'Eliseo D, Ziosi M, Acciarri R. Effects of phacoemulsification time on the corneal endothelium using phacofracture and phaco chop techniques. *J Cataract Refract Surg* 1996; 22:967–969
4. Sippel KC, Pineda R II. Phacoemulsification and thermal wound injury. *Semin Ophthalmol* 2002; 17:102–109
5. Hayashi K, Hayashi H, Nakao F, Hayashi F. Risk factors for corneal endothelial injury during phacoemulsification. *J Cataract Refract Surg* 1996; 22:1079–1084
6. Miyata K, Nagamoto T, Maruoka S, Tanabe T, Nakahara M, Amano S. Efficacy and safety of the soft-shell technique in cases with a hard lens nucleus. *J Cataract Refract Surg* 2002; 28:1546–1550
7. Kim EK, Cristol SM, Geroski DH, McCarey BE, Edelhauser HF. Corneal endothelial damage by air bubbles during phacoemulsification. *Arch Ophthalmol* 1997; 115:81–88; correction, 630
8. Kim EK, Cristol SM, Kang SJ, Edelhauser HF, Kim H-L, Lee JB. Viscoelastic protection from endothelial damage by air bubbles. *J Cataract Refract Surg* 2002; 28:1047–1053
9. Shimmura S, Tsubota K, Oguchi Y, Fukumura D, Suematsu M, Tsuchiya M. Oxiradical-dependent photoemission induced by a phacoemulsification probe. *Invest Ophthalmol Vis Sci* 1992; 33:2904–2907. Available at: <http://www.iovs.org/cgi/reprint/33/10/2904>. Accessed July 14, 2009
10. Holst A, Rolfsen W, Svensson B, Öllinger K, Lundgren B. Formation of free radicals during phacoemulsification. *Curr Eye Res* 1993; 12:359–365
11. Cameron MD, Poyer JF, Aust SD. Identification of free radicals produced during phacoemulsification. *J Cataract Refract Surg* 2001; 27:463–470
12. Takahashi H, Sakamoto A, Takahashi R, Ohmura T, Shimmura S, Ohara K. Free radicals in phacoemulsification and aspiration procedures. *Arch Ophthalmol* 2002; 120:1348–1352. Available at: <http://archophth.ama-assn.org/cgi/reprint/120/10/1348.pdf>. Accessed July 14, 2009
13. Murano N, Ishizaki M, Sato S, Fukuda Y, Takahashi H. Corneal endothelial cell damage by free radicals associated with ultrasound oscillation. *Arch Ophthalmol* 2008; 126:816–821. Available at: <http://archophth.ama-assn.org/cgi/reprint/126/6/816?ck=nck>. Accessed July 14, 2009
14. Suzuki H, Takahashi H, Hori J, Hiraoka M, Igarashi T, Shiwa T. Phacoemulsification associated corneal damage evaluated by corneal volume. *Am J Ophthalmol* 2006; 142:525–528
15. Suzuki H, Oki K, Takahashi K, Shiwa T, Takahashi H. Functional evaluation of corneal endothelium by combined measurement of corneal volume alteration and cell density after phacoemulsification. *J Cataract Refract Surg* 2007; 33:2077–2082
16. Tham CCY, Kwong YYY, Lai JSM, Lam DSC. Effect of a previous acute angle closure attack on the corneal endothelial cell density in chronic angle closure glaucoma patients. *J Glaucoma* 2006; 15:482–485



First author:

Hisaharu Suzuki, MD, PhD

Department of Ophthalmology, Nippon Medical School, Tokyo, Japan